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COMPONENTES E INSTRUMENTOS DE
PROSTODONCIA

EN | USAGE INSTRUCTIONS

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USAGE INSTRUCTIONS

IMPORTANT INFORMATION.

Read this prospectus carefully before using the product.

This document contains information for the use of prosthetic instruments and attachments or components for the KYT® dental implant systems. This documentation has been drafted as a quick guide to be consulted by the final user of the system and does not provide the indications and technical specifications necessary to properly use the products in the KYT® dental implant systems. This is not an alternative nor substitute for specialized training or professional clinical experience.

The KYT® dental implant systems must be used following a strict surgical and prosthetic protocol that does not excessively deviate from any other conventional implant on the market. As with these systems, adequate treatment planning is necessary. However, one must take some differences from other implant systems into account, both in planning treatment and in the surgical and prosthetic protocol in order to avoid errors in the working process and manipulation of the products.

Read the entire surgical and prosthetic protocol, the usage instructions and the maintenance instructions carefully before using any product from the KYT® implant systems. You can consult our web site www.kytimplantesdentales.com or ask your nearest KYT® distributor.

1. Description of the system.

KYT® has developed several dental implant systems in order to provide dental treatments in accordance with the existing scientific literature and current clinical standards.

KYT® develops, manufactures and markets ultimate generation surgical instruments and a wide array of restoration solutions to meet each patient's particular needs.

2. Description of the products.

Instrument

KYT® instruments are manufactured in stainless steel and titanium. The surgical burrs are designed to support a maximum of 45 uses. Inadequate use, lack of cleaning, maintenance or inadequate disinfection may reduce your number of uses. Some of the instruments are identified by laser etching for easy identification and proper use.

Included in the instruments are:

Surgical Burrs - Burr stops - Carriers - Bone Trephines - Bone shapers - Manual or mechanical insertion adapters - Insertion positioners - Ratchets - Manual or mechanical screwdrivers - Meters - Abutment adapters or components, among others.

Attachments

KYT® prosthetic components or attachments are manufactured with titanium, stainless steel and biocompatible plastic materials. Some KYT® components, abutments or screwdrivers are identified using colour codes for easy identification and adequate use depending on each implant series.

Among the prosthetic attachments or components, we can find:

Healing abutments - Impression molds - Straight and angled modifiable abutments - Provisional titanium abutments and biocompatible plastic (Peek) - Basic® transepithelial abutments - Straight and angled Jocker® - O-Ring® abutments - Simpro® abutments - Netcan® abutments - Clinical screws and laboratory screws, among others.

3. Usage indications.

Before performing any treatment with KYT® dental implants, it is essential that you evaluate the patient and the patient's possible contraindications for treatment with dental implants (clinical diagnosis, radiology and study of dental models). Proper treatment with dental implants must first have a proper prosthodontic rehabilitation design that meets the patient's needs and expectations.

Instrument

KYT® surgical instruments are designed as one of the fundamental parts of implant treatment. Using the instruments, you can carry out all of the processes to prepare the bone bed, insert the dental implants and handle the different KYT® prosthetic attachments or components.

Attachments

The KYT® prosthetic attachments or components are attachments or abutments that fix the KYT® dental implants and serve as support for prosthetic rehabilitation. Prosthetic rehabilitation can be single, partial or complete and can be screwed or cemented to the attachments or being removable using implant anchors or retainers.

The biocompatible plastic (Peek) attachments or abutments are for provisional use and are designed for use in the mouth for a period of no longer than 180 days. Removable prosthetics with anchors or retainers must be monitored at least one time per year by the clinical professional. You must verify proper adjustment of the prosthesis in the mouth and retention between the anchors or retainers and the dental prosthesis. Loss of adjustment or poor retention of the prosthetic put prosthodontic rehabilitation and dental implants at risk.

4. Contraindications.

A preoperative medical examination of the patient is necessary to determine if there are contraindications or risk factors for treatment with dental implants. KYT® dental implants should not be used in patients who lack the minimal medical conditions necessary to carry out prosthodontic treatment and rehabilitation with dental implants.

Relative. The clinical professional is responsible for evaluating the potential benefits and risks for the particular treatment: Bone deficiency - Bone quantity/quality - Patient's age - Patient's health - Poor oral hygiene - Localised infections at implant site - Residual roots at implant site - Pregnancy - Stress - Mental factors - Prosthetic heart valves - Bruxism - Peri-dontal pathology - Gum pathology - Tobacco use - Alcohol abuse - Drug Addiction - Addiction in general, among others.

Absolute. The manufacturer does not recommend the treatment in the next cases: Cardiovascular diseases or disorders - Haematological diseases - Maxillary sinusitis - End-stage disease - Osteoporosis or bone disorders - Epilepsy - Neoplasms treated with chemotherapy - Blood dyscrasias - Cancer - Severe or poorly controlled diabetes - Untreated periodontitis

- Infectious diseases - Leukaemia - Immune system alterations - Allergy to titanium - Allergy stainless steel. - Radiation therapy treatment - Corticosteroid treatments - Anticoagulant treatments - Involuntary muscle disorders - Severe liver disorders - Endocrine disorders - Uncontrolled hypertension - Psychiatric and/or psychological disorders.

5. Storage.

KYT® products must be stored in their original packaging in a dry, clean location protected from sunlight and adverse conditions.

6. Packaging and manner supplied.

KYT® instruments or prosthetic attachments or components come in a thermosealed blister pack. The product inside is identified with an identifying label at the back of the blister. In order to do obligatory proper monitoring of the traceability of the products, you must transcribe the patient's identification, the reference number and the lot number of the product being used.

KYT® guarantees that all of its products follow a strict manufacturing, quality control and cleaning process before being packaged. Before using a KYT® implant product, you must check for proper integrity of the packaging and ensure that it is not damaged. If there is a defect or damage, the product must not be used, must be discarded and a replacement product ordered. KYT® will replace any damaged product at no cost so long as it does not show signs of previous handling.

7. Cleaning and sterilization.

The KYT® dental implant systems instruments and prosthetic attachments or components are supplied in UN-sterilized. The instruments and attachments or components that come into contact with the mouth must be cleaned, disinfected and sterilized prior to use.

We recommend: Open the original package and remove the product from its interior, submerge it in an ultrasonic bath with a mixture of detergent and water for at least 10 minutes, remove the product and then deposit it in bags for sterilization in a steam autoclave at 134 °C for at least 5 minutes (UNE ISO 17665-1). In the case of plastic elements, submerge in 20% Chlorhexidine for 10 minutes or sterilize using cold disinfection with chemicals. Afterwards, rinse the products for two minutes with germ and endotoxin-free deionized water.

In reusable products, you must clean, disinfect and sterilize after each use. For more information, consult the specific protocol described in the document CLEANING AND STERILIZATION on our web site www.kytimplantesdentales.com or ask your nearest KYT® distributor.

8. Precautions.

Surgery and Instruments

The surgical protocols describe in detail the precautions and important factors to keep in mind during the surgical process for placing the dental implants. Due to the average dimensions of the products used, among the most important measures you must take special care to ensure that these products are not accidentally ingested or aspirated by the patient. All KYT® instruments for manual use incorporate retention elements for use with dental thread or silk to avoid accidental ingestion.

You must do everything possible to minimize the surgical time and damage to the area where the implants will be placed, paying special attention to bone trauma and asepsis of the surgical field. Proper preparation of the surgical bed is a fundamental factor for successful treatment. This process requires that you follow the surgical burring protocol provided by the manufacturer closely for each type of implant and proper use of the surgical burrs, both in the number of uses and the working speeds. The surgical burrs must be used with constant abundant irrigation to avoid bone overheating.

Pay special attention to not exerting excessive forces or torque during placement of the implant. It is important to control the insertion torque of the implant, both in mechanical insertion and in manual insertion.

Attention. Use of the surgical burrs, Maximum 45 Uses. Surpassing the number of uses of the surgical burrs puts the osseointegration process and therefore the dental treatment at risk. The maximum insertion torque for the dental implants is 50 Ncm. Surpassing the maximum insertion torque for the implants may cause severe damage in the dental implant, its connection, the carrier and the clinical screw.

Prosthodontic Rehabilitation

The prosthetic protocols describe in detail the precautions and important factors to keep in mind during the surgical process for placing prosthetic rehabilitation. The design, type of rehabilitation and placement of the dental prosthetic must be studied and carried out prior to placement of the dental implants. Be careful to not apply excessive force during the tests of the attachments or dental prosthetics. It must always be done manually. Pay special attention when doing the final adjustment of the rehabilitation. You must correctly control the torque applied to the final clinical screws according to the manufacturer's indications whether you do it mechanically or manually.

Adjustment torques for KYT® screws

O-Ring® abutments/ Simpro® / Netcan® / Jocker®- 25 Ncm

Standard Jocker® clinical screws- 20 Ncm

All laboratory screws of KYT systems- 10 Ncm

Attention. Final screws must have the correct indicated torque in order to adequately carry out their purpose. Not applying the correct torque to the clinical screws indicated by the manufacturer, whether by default or due to excess torque, may cause severe damage in the dental implant, its connection, the abutment and the clinical screw. Always confirm passive adjustment of the implant-prosthetic unit. Poor passive adjustment puts the prosthetic rehabilitation in serious risk and may lead to future damage in the dental implant, its connection, the restoration abutment, the clinical screw and/or bone loss or resorption.

9. Warnings.

Treatment planning and placement of KYT® dental implants require specific dental training. The product specifications themselves are not guarantee of proper use. We recommend that clinical users take courses with specific theoretical/practical training in order to understand the proper working techniques and protocols for each product, including the biomechanical, radiographic and prosthetic requirements associated with implant treatment. KYT® periodically gives various courses with which you can update your knowledge on the working techniques and procedures with their implant systems. The end user must be knowledgeable and inform himself about the current state of the technique for each product and its possible applications.

Before using KYT® implant systems or any of their prosthetic attachments or components, it is necessary to understand and be familiar with the corresponding surgical and prosthetic procedures. The patient must meet some sufficient anatomical and psychological conditions to undergo treatment with dental implants. The person responsible for implant treatment must watch for proper treatment planning, guaranteeing the safety of the treatment with a minimal margin of error, respecting the patient's vital mouth structures and general health.

The procedure for using the surgical burrs, saws and other surgical instruments needed to place the implant are specified in the corresponding surgical protocols for each implant. Placement of the implant and prosthodontic planning must adapt to the patient's individual conditions, particularly to proper distribution of forces. Passive adjustment in the prosthodontic rehabilitation, adjustment of occlusal contact points at the opposing maxillary and avoiding the onset of excessive lateral forces must be achieved. Treatment with an insufficient number of implants, improper choice of the size or an inappropriate position to support and transmit the expected loads, may lead to mechanical failure of the implant, the abutment or the prosthetic attachment, the clinical screws due to overload or fatigue and substantial loss of bone around the implant.

Lack of an adequate amount of residual bone, onset of infection or diseases in general and changes in the patient's habits are some potential causes of osseointegration and treatment failure. Lack of bone or soft tissue may lead to unfavorable insertion of the implant and lead to a subsequent deficient aesthetic outcome. Inadequate prosthodontic rehabilitation may cause failure of the entire rehabilitation treatment.

Reuse of single-use products carries the risk of possible deterioration of their characteristics. The geometry of the product may have suffered damage from previous use or due to inadequate use. This involves the risk of surgical or prosthetic failure and/or damage to the patient's health.

Electrosurgery is not indicated due to the conductivity of dental implants.

10. Adverse or side effects.

Current implant technology can have several adverse or undesired effects that are documented in the specialized scientific literature published in journals or specialized books on dentistry. However, the most relevant are:

Postoperative malaise- Local pain- Local inflammation- Hematomas- Local or systemic infections- Difficulty speaking- Bone loss and/or fractures- Loss of the implant- Damage to adjacent teeth - Fractures of implants and prosthetic components - Instrument fracture- Damage to the tooth nerve- Damage to the maxillary sinuses, among others.

11. Information on responsibility, safety and guarantee.

Our usage and application indications for our products are based on the internationally published literature, current clinical standards and our clinical experience with our dental implant system, so they must be understood as general information on indications. Handling and use of the KYT® dental implant systems, being out of the control of KYT®, are the sole responsibility of the end user of the products. KYT Implantes Dentales, SL and the distributor of the KYT® implant systems decline all civil responsibility for possible damage or injury caused by poor handling of the product or any other situation not considered in the protocol for proper use of its products.

The end user must ensure that the product used is adequate for the planned purposes and procedures. Neither these usage instructions, the working protocols or manipulation of the products release the user from this obligation. The KYT® dental implant systems must be used, handled and applied by professionals with the necessary qualifications required according to current legislation of each country. Total or partial use, handling or application of any phase of the KYT® dental implant systems by non-qualified personnel or by personnel who do not have the necessary training will automatically void any type of warranty and may cause serious damage to the patient's health.

The KYT® dental implant systems products are part of their own system, with their own design characteristics and working protocols that include dental implants, prosthodontic attachments or components and surgical or prosthetic instruments. The use of original KYT® products in combination with inappropriate elements or components or products from other manufacturers may lead to failure of the rehabilitation treatment, cause tissue damage, lead to damage to bone structures, inadequate aesthetic outcomes and severe damage to the patient's health.

For this reason, ONLY original KYT® products should be used. The clinical professional is responsible for ensuring the use of original KYT® products and that they are used in compliance with the usage instructions and the corresponding handling protocols throughout the implant treatment process. The use of non-genuine KYT® components, instruments or any other part or element, whether used alone or in combination with any genuine KYT® products, will void any full or partial warranty of all genuine KYT® implants products.

In order to receive the guaranty for the original KYT® products, all of the products used in the implant treatment must be original KYT® products and must meet the planned usage requirements for the product, the usage instructions and the original KYT® handling protocols. Both the person responsible for the dental treatment and the patient must provide all information required by KYT® for study. All original products that have been used to treat patients with risk factors or contraindications expressed in this document are excluded from the total or partial guarantee.

Warning. Not all KYT® products are for sale in every country.

For more information, consult the specific surgical or prosthesis protocols for each product at our web site www.kytimplantesdentales.com or ask your nearest KYT® distributor.

Ktx®, Kti® and all their logos are registered trademarks of KYT®.

The KYT® implant systems are a registered trademark, as are some of the products mentioned or not mentioned in this catalogue. KYT Medical, SL reserves the right to modify, change, eliminate or update any of the products, prices or technical specifications referenced in this document or in any of its catalogues without prior notification. All rights reserved. Any complete or partial reproduction of this documentation, in any medium or format, is forbidden without the corresponding written authorization from KYT Implantes Dentales, SL

12. Magnetic resonance.

The KYT® dental implant systems products have not been evaluated for safety and compatibility nor have they been tested for heating and migration in the magnetic resonance environment.

13. Measuring function.

Devices with a measuring function provide pinpoint accuracy.

When you work with burrs you must take into account the length of the burrs tip (0.5 mm) as it is not included in the measurements of the different laser brands.

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