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COMPONENTES E INSTRUMENTOS
DE PROSTODONCIA

EN | INSTRUCTIONS FOR USE

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USAGE INSTRUCTIONS

IMPORTANT INFORMATION

Read this prospectus carefully before using the product.

This document contains information for the use of prosthetic instruments and abutments or components for the KYT® dental implant systems. This documentation has been drafted as a quick guide to be consulted by the final user of the system and does not provide the indications and technical specifications necessary to properly use the products in the KYT® dental implant systems. This is not an alternative nor substitute for specialized training or professional clinical experience.

The KYT® dental implant systems must be used following a strict surgical and prosthetic protocol that does not excessively deviate from any other conventional implant on the market. As with these systems, adequate treatment planning is necessary. However, one must take some differences from other implant systems into account, both in planning treatment and in the surgical and prosthetic protocol in order to avoid errors in the working process and manipulation of the products.

Read the entire surgical and prosthetic protocol, the usage instructions and the maintenance instructions carefully before using any product from the KYT® implant systems. You can consult our web site www.kytimplantesdentales.com or ask your nearest KYT® distributor.

1. Description of the system

KYT® has developed several dental implant systems in order to provide dental treatments in accordance with the existing scientific literature and current clinical standards.

KYT® develops, manufactures and markets ultimate generation surgical instruments and a wide array of restoration solutions to meet each patient's particular needs.

2. Description of the products

Instrument

KYT® instruments are manufactured in stainless steel and titanium. Some of the instruments are identified by laser etching for easy identification and proper use.

Included in the instruments are:

Surgical Burrs - Burr stops - Carriers - Bone Trephines - Bone shapers - Manual or mechanical insertion adapters - Insertion positioners - Ratchets - Manual or mechanical screwdrivers - Probes - Abutment adapters or components, among others.

Abutments

KYT® prosthetic components or abutments are manufactured with titanium, stainless steel and biocompatible plastic materials. Some KYT® components, abutments or screwdrivers are identified using colour codes for easy identification and adequate use depending on each implant series.

Among the prosthetic abutments or components, we can find:

Healing abutments - Impression molds - Straight and angled modifiable abutments - Provisional titanium abutments and biocompatible plastic (Peek) - Basic transepithelial abutments - Straight and angled Jocker - O-Ring abutments - Simpro abutments - Netcan abutments - Clinical screws and laboratory screws, among others.

3. Usage indications

Before performing any treatment with KYT® dental implants, it is essential that you evaluate the patient and the patient's possible contraindications for treatment with dental implants (clinical diagnosis, radiology and study of dental models). Proper treatment with dental implants must first have a proper prosthodontic rehabilitation design that meets the patient's needs and expectations.

Instrument

KYT® surgical instruments are designed as one of the fundamental parts of implant treatment. Using the instruments, you can carry out all of the processes to prepare the bone bed, insert the dental implants and handle the different KYT® prosthetic abutments or components.

The reusable surgical burs and trephines are designed to be used up to 45 times under optimum conditions; whilst the circular scalpels are designed to withstand a maximum of 25. Improper use, lack of cleaning and maintenance or inadequate disinfection can reduce its number of uses.

The KYT® instruments are re-usable.

Abutments

KYT® abutments and prosthetic components are designed to be connected to dental implants and are intended to support prosthetic rehabilitation. Prosthetic rehabilitation can be single, partial or complete and can be screwed or cemented to the abutments or being removable using implant anchors or retainers.

The biocompatible plastic (Peek) abutments or abutments are for provisional use and are designed for use in the mouth for a period of no longer than 180 days. Removable prosthetics with anchors or retainers must be monitored at least one time per year by the clinical professional. You must verify proper adjustment of

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the prosthesis in the mouth and retention between the anchors or retainers and the dental prosthesis. Loss of adjustment or poor retention of the prosthetic put prosthodontic rehabilitation and dental implants at risk.

The KYT® abutments are single-use products.

4. Intended patients

Instruments

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment in the mandible and maxilla, and are intended to restore chewing function, as well as providing an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

The target population age, approximately, can be from 6 years old with orthodontic treatment.

Abutments

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment in the mandible and maxilla, and are intended to restore chewing function, as well as providing an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

These patients must have completed their growth. The age guideline is older than 18 years, this can be confirmed after cessation of growth with a wrist X-ray. Patients at an early age patients will be specifically informed of the changes they may need to make to the prostheses during their lifetime with a with the consequent aesthetic and economic impact due to functional ankylosis of the implant in living tissue.

5. Contraindications

A preoperative medical examination of the patient is necessary to determine if there are contraindications or risk factors for treatment with dental implants. KYT® dental implants should not be used in patients who lack the minimal medical conditions necessary to carry out prosthodontic treatment and rehabilitation with dental implants.

Relative

The clinical professional is responsible for evaluating the potential benefits and risks for the particular treatment: Bone deficiency - Bone quantity/quality - Patient’s age - Patient’s health - Poor oral hygiene - Localised infections at implant site - Residual roots at implant site - Pregnancy - Stress - Mental factors - Prosthetic heart valves - Bruxism - Peridental pathology - Gum pathology - Tobacco use - Alcohol abuse - Drug Addiction - Addiction in general, among others.

Absolute

The manufacturer does not recommend the treatment in the next cases: Cardiovascular diseases or disorders - Haematological diseases - Maxillary sinusitis - End-stage disease - Osteoporosis or bone disorders - Epilepsy - Neoplasms treated with chemotherapy - Blood dyscrasias - Cancer - Severe or poorly controlled diabetes- Untreated periodontitis - Infectious diseases - Leukaemia - Immune system alterations - Allergy to titanium - Allergy stainless steel. - Radiation therapy treatment - Corticosteroid treatments - Anticoagulant treatments - Involuntary muscle disorders - Severe liver disorders - Endocrine disorders - Uncontrolled hypertension - Psychiatric and/or psychological disorders.

6. Storage

KYT® products must be stored in their original packaging in a dry, clean location protected from sunlight and adverse conditions.

7. Packaging and manner supplied

KYT® instruments or prosthetic abutments or components come in a thermosealed blister pack. The product inside is identified with an identifying label at the back of the blister. In order to do obligatory proper monitoring of the traceability of the products, you must transcribe the patient’s identification, the reference number and the lot number of the product being used.

Included in the abutments, the removable labels serve as a traceability control of the abutment. This is achieved by sticking one label to the implant card for the patient, and the rest for the clinic’s record.

KYT® guarantees that all of its products follow a strict manufacturing, quality control and cleaning process before being packaged. Before using a KYT® implant product, you must check for proper integrity of the packaging and ensure that it is not damaged. If there is a defect or damage, the product must not be used, must be discarded and a replacement product ordered. KYT® will replace any damaged product at no cost so long as it does not show signs of previous handling.

8. Cleaning and sterilization

The KYT® dental implant systems instruments and prosthetic abutments or components are supplied in UN-sterilized. The instruments and abutments or components that come into contact with the mouth must be cleaned, disinfected and sterilized prior to use.

9. Precautions

Surgery and Instruments

The surgical protocols describe in detail the precautions and important factors to keep in mind during the surgical process for placing the dental implants. Due to the average dimensions of the products used, among the most important measures you

must take special care to ensure that these products are not accidentally ingested or aspirated by the patient. All KYT® instruments for manual use incorporate retention elements for use with dental thread or silk to avoid accidental ingestion.

You must do everything possible to minimize the surgical time and damage to the area where the implants will be placed, paying special attention to bone trauma and asepsis of the surgical field. Proper preparation of the surgical bed is a fundamental factor for successful treatment. This process requires that you follow the surgical burring protocol provided by the manufacturer closely for each type of implant and proper use of the surgical burs, both in the number of uses and the working speeds. The surgical burs must be used with constant abundant irrigation to avoid bone overheating.

Pay special attention to not exerting excessive forces or torque during placement of the implant. It is important to control the insertion torque of the implant, both in mechanical insertion and in manual insertion.

Attention

Use of the surgical burs, Maximum 45 Uses. Surpassing the number of uses of the surgical burs puts the osseointegration process and therefore the dental treatment at risk.

The maximum insertion torque for the dental implants is 50 Ncm. Surpassing the maximum insertion torque for the implants may cause severe damage in the dental implant, its connection, the carrier and the clinical screw.

When using surgical drills, you must take into account the length of the drill tip (+0.5 mm), as IT IS NOT INCLUDED in the length of the implant to be drilled.

Prosthodontic Rehabilitation

The prosthetic protocols describe in detail the precautions and important factors to keep in mind during the surgical process for placing prosthetic rehabilitation. The design, type of rehabilitation and placement of the dental prosthetic must be studied and carried out prior to placement of the dental implants. Be careful to not apply excessive force during the tests of the abutments or dental prosthetics. It must always be done manually. Pay special attention when doing the final adjustment of the rehabilitation. You must correctly control the torque applied to the final clinical screws according to the manufacturer’s indications whether you do it mechanically or manually.

KYT® torque table for abutments:

Element/Abutment	Instrument/Tool	Torque
Cover screws/healing abutments	Hex. screwdriver 1.25mm	Manual
Print transfer screws	Hex. screwdriver 1.25mm	Manual
Laboratory screws	Hex. screwdriver 1.25mm	10 Ncm
Clinical screws	Hex. screwdriver 1.25mm	30 Ncm
Scanbody + screw	Hex. screwdriver 1.25mm	Manual
Jocker abutments	Key KMA20	30 Ncm
Transocclusal screws on Jocker	Hex. screwdriver 1.25mm	20 Ncm
Simpro abutments	Key KLOSD10 (square 2.11mm)	30 Ncm

Some instruments are designed for use with a contra-angle. In order to guarantee a correct connection of both products, please see the contra-angle manufacturer’s instructions for use.

Attention

Final screws must have the correct indicated torque in order to adequately carry out their purpose. Not applying the correct torque to the clinical screws indicated by the manufacturer, whether by default or due to excess torque, may cause severe damage in the dental implant, its connection, the abutment and the clinical screw. Always confirm passive adjustment of the implant-prosthetic unit. Poor passive adjustment puts the prosthetic rehabilitation in serious risk and may lead to future damage in the dental implant, its connection, the restoration abutment, the clinical screw and/or bone loss or resorption.

10. Warnings

Treatment planning and placement of KYT® dental implants require specific dental training. The product specifications themselves are not guarantee of proper use. We recommend that clinical users take courses with specific theoretical/practical training in order to understand the proper working techniques and protocols for each product, including the biomechanical, radiographic and prosthetic requirements associated with implant treatment. KYT® periodically gives various courses with which you can update your knowledge on the working techniques and procedures with their implant systems. The end user must be knowledgeable and inform himself about the current state of the technique for each product and its possible applications.

Before using KYT® implant systems or any of their prosthetic abutments or components, it is necessary to understand and be familiar with the corresponding surgical and prosthetic procedures. The patient must meet some sufficient anatomical and psychological conditions to undergo treatment with dental implants. The person responsible for implant treatment must watch for proper treatment planning, guaranteeing the safety of the treatment with a minimal margin of error, respecting the patient’s vital mouth structures and general health.

The procedure for using the surgical burs, saws and other surgical instruments needed to place the implant are specified in the corresponding surgical protocols for each implant. Placement of the implant and prosthodontic planning must adapt

to the patient's individual conditions, particularly to proper distribution of forces. Passive adjustment in the prosthodontic rehabilitation, adjustment of occlusal contact points at the opposing maxillary and avoiding the onset of excessive lateral forces must be achieved. Treatment with an insufficient number of implants, improper choice of the size or an inappropriate position to support and transmit the expected loads, may lead to mechanical failure of the implant, the abutment or the prosthetic abutments, the clinical screws due to overload or fatigue and substantial loss of bone around the implant.

Lack of an adequate amount of residual bone, onset of infection or diseases in general and changes in the patient's habits are some potential causes of osseointegration and treatment failure. Lack of bone or soft tissue may lead to unfavorable insertion of the implant and lead to a subsequent deficient aesthetic outcome. Inadequate prosthodontic rehabilitation may cause failure of the entire rehabilitation treatment.

Reuse of single-use products carries the risk of possible deterioration of their characteristics. The geometry of the product may have suffered damage from previous use or due to inadequate use. This involves the risk of surgical or prosthetic failure and/or damage to the patient's health.

Electrosurgery is not indicated due to the conductivity of dental implants.

11. Adverse or side effects

Current implant technology can have several adverse or undesired effects that are documented in the specialized scientific literature published in journals or specialized books on dentistry.

However, the most relevant are:

Postoperative malaise - Local pain - Local inflammation - Hematomas - Local or systemic infections - Difficulty speaking - Bone loss and/or fractures - Loss of the implant - Damage to adjacent teeth - Fractures of implants and prosthetic components - Instrument fracture - Damage to the tooth nerve - Damage to the maxillary sinuses, among others.

12. Information on responsibility, safety and guarantee

Our usage and application indications for our products are based on the internationally published literature, current clinical standards and our clinical experience with our dental implant system, so they must be understood as general information on indications. Handling and use of the KYT® dental implant systems, being out of the control of KYT®, are the sole responsibility of the end user of the products. KYT Implantas Dentales, SL and the distributor of the KYT® implant systems decline all civil responsibility for possible damage or injury caused by poor handling of the product or any other situation not considered in the protocol for proper use of its products.

The end user must ensure that the product used is adequate for the planned purposes and procedures. Neither these usage instructions, the working protocols or manipulation of the products release the user from this obligation. The KYT® dental implant systems must be used, handled and applied by professionals with the necessary qualifications required according to current legislation of each country. Total or partial use, handling or application of any phase of the KYT® dental implant systems by non-qualified personnel or by personnel who do not have the necessary training will automatically void any type of warranty and may cause serious damage to the patient's health.

The KYT® dental implant systems products are part of their own system, with their own design characteristics and working protocols that include dental implants, prosthodontic abutments or components and surgical or prosthetic instruments. The use of original KYT® products in combination with inappropriate elements or components or products from other manufacturers may lead to failure of the rehabilitation treatment, cause tissue damage, lead to damage to bone structures, inadequate aesthetic outcomes and severe damage to the patient's health.

For this reason, ONLY original KYT® products should be used. The clinical professional is responsible for ensuring the use of original KYT® products and that they are used in compliance with the usage instructions and the corresponding handling protocols throughout the implant treatment process. The use of non-genuine KYT® components, instruments or any other part or element, whether used alone or in combination with any genuine KYT® products, will void any full or partial warranty of all genuine KYT® implants products.

In order to receive the guaranty for the original KYT® products, all of the products used in the implant treatment must be original KYT® products and must meet the planned usage requirements for the product, the usage instructions and the original KYT® handling protocols. Both the person responsible for the dental treatment and the patient must provide all information required by KYT® for study. All original products that have been used to treat patients with risk factors or contraindications expressed in this document are excluded from the total or partial guarantee.

Any problem or incident with regard to the device should be reported to the manufacturer, to which the user should also send the affected product. In the event of a serious incident, the user should also submit a report to the manufacturer, as well as the relevant competent authority, according to the requirements of local legislation.

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

Warning

Not all KYT® products are for sale in every country.

For more information, consult the specific surgical or prosthesis protocols for each product at our web site www.kyimplantesdentales.com or ask your nearest KYT® distributor.

Ktx®, Kti® and all their logos are registered trademarks of KYT®.

The KYT® implant systems are a registered trademark, as are some of the products mentioned or not mentioned in this catalogue. KYT Medical, SL reserves the right to modify, change, eliminate or update any of the products, prices or technical specifications referenced in this document or in any of its catalogues without prior notification. All rights reserved. Any complete or partial reproduction of this documentation, in any medium or format, is forbidden without the corresponding written authorization from KYT Implantas Dentales, SL.

13. Magnetic resonance

The KYT® dental implant systems products have not been evaluated for safety and compatibility nor have they been tested for heating and migration in the magnetic resonance environment.

14. Removal of the material

The dental implants, abutments and kit, as well as each product and consumable used during surgery for placing dental implants and which are disposed of may pose a health hazard to people handling them. For this reason, you are advised to check and comply with current legislation on biohazardous waste.

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EN | CLEANING, DISINFECTION AND STERILISATION

Cleaning, disinfection and sterilisation

The protocols described in this section must only be carried out by personnel qualified to clean, disinfect and sterilise the dental materials specified here in.

CLEANING AND DISINFECTION INSTRUCTIONS

Applicable for instruments, surgical and prosthetic boxes.

Disassembly

1. Dismount* the appropriate instruments, for example manual ratchets, drills or drill stops.
2. Remove the various components from the surgical or prosthetic box for correct cleaning.

Cleaning and disinfection

For disinfecting instruments and surgical boxes:

1. Submerge the instruments in a detergent/disinfectant solution** suitable for dental instruments to help eliminate any adhered biological residues. If an ultrasound bath is available***, confirm that the detergent/disinfectant solution is indicated for use with this type of equipment.
2. Manually remove any biological residues with a non-metallic brush and pH-neutral detergent.
3. Rinse with copious water.
4. When cleaning the surgical and prosthetic boxes, always use a pH-neutral detergent and non-abrasive utensils to avoid damaging the surface of the boxes.
5. Dry the materials with disposable cellulose, lint-free clothes or compressed air.

For disinfecting plastic caps and spacers:

1. Submerge in a neat benzalkonium chloride solution for 10 minutes.
2. Rinse with distilled water.
3. Dry the caps and spacer before use.

Inspection

1. Check that the instruments are perfectly clean; if not, repeat the cleaning and disinfection steps.
2. Discard any instruments with imperfections and replace them before the next procedure.

3. Check that the instruments and the surgical and prosthetic boxes are perfectly dry before reassembling the parts and proceeding to their sterilisation.

* See assembly and disassembly instructions.

** Follow the instructions from the disinfectant's manufacturer to determine the correct concentrations and times.

*** Follow the instructions from the ultrasound bath's manufacturer to determine the correct temperature, concentration and times.

STERILISATION INSTRUCTIONS FOR STEAM AUTOCLAVE

Applicable to orthodontic implants, abutments, instruments, surgical and prosthetic boxes.

1. Introduce each material separately in individual sterilisation bags, then seal the bags. For joint sterilisation, place the instruments in their surgical box, introduce the box into a sterilisation bag and seal the bag.
2. Place the bags to be sterilised in the autoclave.
3. Sterilise in a steam autoclave at 134°C/273°F (max. 137°C/276°F) for 4 min (minimum) and at 2 atm. Torque wrenches must be sterilised in 3 vacuum cycles at 132°C/270°F for a minimum of ≥ 4 minutes and vacuum-dried for a minimum of 20 minutes.

For the United States only: The validated and recommended sterilisation cycle for the US must be performed in a steam autoclave at 132°C/270°F for at least 15 min and with the drying time of at least 15 - 30 min.

Important:

- Make sure the drying stage is allowed to run to completion, otherwise the products may be damp.
- Check the sterilisation equipment if the materials or sterilisation bags are damp at the end of the sterilisation cycle.
- Perform the necessary maintenance actions on the autoclave according to the established periodicity and following the manufacturer's instructions.

STORAGE OF KYT® PRODUCTS

- Store the products in their original packaging and in a clean, dry location until they are used.
- After sterilisation, keep the products in the sealed sterilisation bags and in a clean, dry location.
- Never exceed the use by date indicated by the manufacturer of the sterilisation bags.
- Always follow the indications of the manufacturer of the sterilisation bags.

GENERAL RECOMMENDATIONS

- Never use damaged or dirty material; never reuse single-use products. The user is responsible for following the instructions described in this document correctly.
- The attention to piercing or sharp elements. Gloves should be worn when cleaning the materials to avoid accidents during handling.
- Follow the safety instructions indicated by the manufacturer of the disinfectant agent.
- The product's sterility cannot be guaranteed if the sterilisation bag is open, damaged or damp.
- Respect all stages of the sterilisation process. If the materials or sterilisation bags contain traces of water or moisture, check the autoclave and repeat the sterilisation.
- Orthodontic abutments and implants are supplied UNSTERILISED and must always be sterilised before use.
- Instruments and surgical and prosthetic boxes are supplied UNSTERILISED and must always be sterilised before use and cleaned and disinfected after use.
- The sterilisation, cleaning and disinfection processes gradually deteriorate the instruments. Inspect the instruments thoroughly to detect any signs of deterioration.
- Avoid contact between products made from different materials (steel, titanium, etc.) during the cleaning, disinfection and sterilisation processes.
- Ziocom Medical SL recommends these instructions are implemented for the correct maintenance and safety of their products; accordingly, the company refuses any liability for any damage to the products that could arise if the user applies alternative cleaning, disinfection and sterilisation procedures.

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